

# **PERICARDIAL ASPIRATION KIT**

Aspiration set to puncture the pericardium by a Pigtail catheter

The catheter has additionally to the distal port 8 side holes for secretion suction, as well as a device for stretching the pigtail tip to make the inserting process easier.



REF	Description
509 045	Pigtail catheter 5 F, 45 cm with three-way-stopcock
	18 G Seldinger cannula 7 cm
	Guide wire J-tip 90 cm x 0.032" with insertion aid
	Suction bag 2.0 ltr. with LL-adaptor and 90 cm suction tube
509 095	Pigtail catheter 5 F, 90 cm
	18 G Seldinger cannula 7 cm + 19 cm
	Dilator 6F - 10 cm
	Guide wire J-tip 150 cm x 0.032" with insertion aid
Box: 5 pcs.	
010 0 <b>45</b>	Collection bag 1.5 ltr. with 90 cm tube
Box: 20 pcs.	

Other lengths and sizes on enquiry.



## **Endovenous Laser Therapy**

Varicose veins pose a problem to approx. 15-30% of the adult population, increasing with age. In milde cases patients are bothered by the unaesthetic appearance of varicose veins. Severe cases are associated with pain, inflammation and in rare cases open wounds. Following needle positioning an introducer is easily inserted into the vein by Seldinger technique. After insertion of a laser fibre the treatment can begin. With the laser fibre in correct position of the saphenofemoral junction (determined by ultrasound) you are ready to deliver the first laser pulse. The laser energy will result in permanent thermal damage to the vein wall. The laser treatment is performed along the entire length of the varicose vein and will result in a complete and permanent occlusion of the vein. After a treatment of 30 – 45 minutes the patient can return to normal activity level.



### Advantages of the method:

- minimal risk of side effect
- better aesthetic results compared to other methods
- no post-op downtime for patients
- · fast out-patient procedure
- only local anaesthesia
- · cost-effective solution

We offer special introducers for bare fibers and radial fibers.

The kit contains the catheter, a puncture needle, dilator und guide wire.

REF	Description	Catheter-Specifications	ID (mm)						
Introducers for bare fibers.									
331 704	Catheter Kit	6F – 70 cm/with extension line, blue soft tip	0,92						
331 004	Catheter Kit	6F – 100 cm/with extension line, blue soft tip	0,92						
330 7055	Catheter Kit	5F – 70 cm; stiff, no extension	1,05						
330 005S	Catheter Kit	5F – 100 cm; stiff, no extension	1,05						
Box: 25 pcs.									

010 015 Lock-Adapter

Box: 50 pcs.

#### Introducers for radial fibers.

300 806 Introducer Kit 6 F 300 807 Introducer Kit 7 F Box: 10 pcs.

Further versions on enquiry



### Index

REF	Page	REF	Page	REF	Page	REF	Page
010 045	24/26	323	<b>9</b> /10	503 002	2 <b>2</b>	503 201	2 <b>2</b>
010 015	27	330	27	503 003	2 <b>2</b>	503 211	2 <b>2</b>
		331	8	503 010	24	503 401	2 <b>2</b>
011	1 <b>3</b> /2 <b>4</b>	331 004	27	503 011	2 <b>2</b>	503 501	2 <b>2</b>
022	16	331 704	27	503 012	2 <b>2</b>	504	24
052 225	16	332	<b>9/</b> 10	503 014	2 <b>2</b>		
100	18	334	8			509	26
103	19	401	13	503 019	23	611 010	13
		402	13	503 021	23	621	13
		440	10	503 022	23		
120	18	451	13	503 025	23	701	14
250	1 <b>0</b>	452	13	503 050	23		
252	10	470	12	503 060	2 <b>3</b>	801	13
300	12	471	12	503 061	23		
302	6	472	12	503 111	2 <b>3</b>		
305	6	473	12	503 119	23		
		503 001	2 <b>2</b>	503 180	24		

## **Information**

All products are delivered **sterile** (ethylene oxide), if not otherwise described.

The products are for **single use only**, if not otherwise described.

Most of intra-products are **latex-free**. But for some products of our range it's not possible to do without.

This is for the following products:

Flow-directed balloon-tipped catheter (Intrathermodin)

Pneumovent® valve

Catheters for vessel surgery (Embolectomies)

Balloon from Latex

**→** Inside lying valve from Latex

Balloon from Latex

