

The image displays several medical catheters and components against a light blue background. At the top, a long, thin, blue catheter is curved. Below it, a red catheter with a red handle and a black catheter with a black handle are visible. In the lower half, there is a large blue T-shaped handle connected to a long, thin, silver catheter. To the right, a smaller blue catheter is partially visible. The text 'Elektrophysiologie' is centered in the middle of the image.

# Elektrophysiologie

# Intralektrode Bipolar

Bipolare Semi-Einschwemm-Elektroden in 4, 5 und 6 French zur temporären Stimulation des Herzens besonders bei Notfällen wie z.B. bradycarde Herzrhythmus- und Reizleitungsstörungen und zur Ableitung des intrakardialen EKG. Die Intralektrode kann schnell, leicht und sicher positioniert werden und zeichnet sich durch gute Vorschubeigenschaften, gute Torsionsübertragung und gute Andruckelastizität aus. Dadurch erreicht der Anwender schnell die gewünschte Signalübertragung. Die Intralektrode lässt sich leicht und stabil vorformen.



## Inhalt

### Intralektrode bipolar

- 1 Stimulationselektrode, 110 cm
- 1 Einführnadel

### Intralektrode mit Ballon

- 1 Stimulationselektrode mit Ballon, 110 cm
- 1 Spritze für den Ballon

REF	French Größe	Spitze
<b>Intralektrode bipolar</b>		
022 244M	4	C-geformt
022 245M	5	C-geformt
022 246M	6	C-geformt
<b>VPE: 5 Stück</b>		
<b>Intralektrode mit Ballon</b>		
052 225	5	Gerade
<b>VPE: 5 Stück</b>		

Elektroden in gerader Ausführung oder mit J-Spitze auf Anfrage.